

DaVinci

Poultry Vaccinator

0.5 Vial Holder Set - Utility

With Plastic Reservoir



Catlog No.: 71-02-22-30

Also Available : 1ml, 2ml, 5ml

0.5 DaVinci Glass Barrel
Catalog No.: 75-02-00-50



0.5 DaVinci Spare Kit
Catalog No.: 72-02-00-45

- Spring, Ball,
- Vial Washer,
- Barrel Washer,
- Intake Pusher &
- Intake Washer



Features

- Plunger tightness without 'O' Ring.
- Barrel made of Borosilicate glass
- Clear visibility of Dose
- Plunger made of Stainless Steel
- Comfortable grip
- Simple and natural movements
- Robust construction
- High quality construction

Advantages

- Fast and reliable volume selection
- Easy changing of dose without changing any piston
- No need for servicing